



## DeVilbiss VacuAide® Compact Suction Unit

### Unmatched safety and performance.

For nearly 50 years, DeVilbiss Healthcare has set the standard of excellence in the homecare suction device category. Now these safety and performance advantages are built into the compact DeVilbiss VacuAide that is designed to fit the active lifestyles of today's customers.

In fact, the VacuAide Compact is the smallest and strongest suction unit of its kind, while featuring a convenient carrying bag that provides discrete suction on the go.



DeVilbiss VacuAide  
Compact Suction Unit

### Performance Features

- Vacuum adjustment allows for 50–550 mm Hg and a free flow of 27 liters/minute
- Meets American Association for Respiratory Care (AARC) guidelines for neonatal, infant, child and adult homecare suctioning
- Enhanced battery life for up to one-hour operation
- Meets ISO10079-1:1999 high vacuum, high flow classification

### Design Features

- Compact package weighs only 3.37 lb
- Easy-to-read LED bar graph gauge helps caregivers provide precise measurements of vacuum – especially important for pediatric patients
- Built-in bacteria filter is located in the unit's bottle assembly; unique float valve prevents backflow into the unit
- Travel-friendly design meets RTCA/DO-160D Aircraft Standard for Environmental Conditions and Test Procedures for Airborne Equipment (won't interfere with a plane's navigation system)

## Design Features

Manufactured and tested to ISO10079-1:1999 Standard for Suction Equipment, achieving outstanding performance in the following tests:

- Electromagnetic Compatibility Test – VacuAide does not interfere with or is not susceptible to interference from other equipment such as ventilators
- Drop Test and Shock Vibration Test – VacuAide is far less likely to be damaged when dropped or shipped
- Excessive Temperature Test – VacuAide effectively works in a wide range of environmental temperatures
- Collection Canister and Tubing Collapsibility Test – The DeVilbiss reusable, autoclavable canister provides reliable, consistent performance



DeVilbiss VacuAide Compact Suction – Model 7310PR-D

## Accessories/Parts

7310P-603	Reusable Collection Canister (bottle, lid, filter, elbow)
7305D-608	Bacteria Filters, 12/pack
7305P-613	AC to DC Adapter/Charger
6305D-611	6' Patient Tubing
7304D-619	12V DC Power Cord
7310P-601	Battery Pack
7310P-602	Battery Door
7310P-606	Carrying Bag

## Technical Specifications

Built to meet the requirements of HCPCS codes:  
**E0600** (Suction Unit), **A7001** (Reusable Collection Bottle)

<b>Product Dimensions</b>	7310PR-D 7.25" h x 7.25" w x 6.75" d	<b>Collection Bottle Capacity</b>	7310PR-D Reusable Autoclavable 725 ml
<b>Product Weight</b>	7310PR-D 3.37 lb	<b>Electrical Requirements</b>	100-240 VAC, 50/60 Hz, 0.75 A max.; 12 V DC, 33 W max.
<b>Shipping Dimensions</b>	7310PR-D 9.25" h x 9.5" w x 12" d	<b>Approvals</b>	IEC 601-1, CAN/CSA-C22.2 No. 601.1-M90, UL 2601-1, CE EN 60601-1-2, ISO10079-1:1999
<b>Shipping Weight</b>	7310PR-D 7.75 lb	<b>Power Consumption</b>	33 watts maximum
<b>Vacuum Range</b>	50-550 mm Hg; meets American Association for Respiratory Care guidelines for neonatal, infant, and adult homecare suctioning	<b>Warranty</b>	2-year warranty on VacuAide unit 6-month warranty on battery
<b>Flow Range</b>	27 LPM free flow, unrestricted		

Conditions may vary based on altitude above sea level, barometric changes and system configuration. Features and specifications are subject to change without notice.



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)