

## EasyFit® SilkGel Nasal Mask

### Lightweight Mask with Superior SilkGel Comfort

The EasyFit SilkGel Nasal Mask has all the same great features of the original EasyFit Nasal Mask, but the mask and forehead cushion have an ultra-soft coating. The smooth surface is soft and silky for exceptional comfort. In addition, the EasyFit SilkGel **features redesigned headgear for a more secure fit and ease of use.**

The EasyFit Nasal Mask was designed by utilizing computer analyses of facial characteristics from 50,000 individuals globally. Identifying consistencies between the profiles attributed to the development of the **unique, scallop-shaped cushion.**

The ball-and-socket joint at the tubing connector **accommodates even the most active sleeper** and contributes to the mask's overall comfort. The unique whisper-quiet EasyFit exhalation system is designed in such a way that neither the patient nor sleeping partner is disturbed by the flow of exhaled air.



EasyFit SilkGel  
Nasal Mask

### Performance Features

- **Exceptionally quiet operation**
- Extremely easy to clean – perspiration and facial oils will not adhere to the smooth cushion and forehead surface
- New mask frame and head gear design
- Extremely durable and intuitive to size
- Unique ball-and-socket joint allows for unsurpassed freedom of movement and a sure and stable mask fit

## Design Features

- Ultra-soft replaceable SilkGel cushions and forehead pads
- Color-coded Velcro fasteners make the headgear even easier to manage
- Tube anchoring strap on the headgear as optional fastening option
- Available in Small, Medium, and Large sizes
- Adjustable forehead support with wider than average forehead cushions for **improved customization**
- Nose flange to prevent mask leaks near the eyes
- Quiet exhalation system **allows patient and bed partner to sleep soundly**
- Multi-patient use approved
- Latex-free



Unique crossing headgear adds stability and minimizes leaks due to active sleeping.



Ultra-soft finish on mask cushion and forehead pad for maximum comfort.

## Reimbursement Information

| Product Description | Code Description   | HCPCS Coding | Frequency      | Allowable Range                  |
|---------------------|--|--------------|----------------|----------------------------------|
| Mask Frame          | Nasal interface (mask or cannula type) used with positive airway pressure device, with or without head strap | A7034        | 1 per 3 months | \$92.57 – \$108.90 initial setup |
| Headgear            | Headgear used with PAP device, each  | A7035        | 1 per 6 months | \$31.27 – \$36.79                |
| Cushion             | Cushion for use on nasal mask interface, replacement only, each  | A7032        | 2 per 1 month  | \$31.89 – \$37.52                |

The noted allowable rates are based off of the January 2012 HCPCS Fee Schedule. All information presented is for reference purposes only.

## Ordering Information

### PRODUCT DESCRIPTION

### ITEM NUMBER

#### Complete System

#### Frame, Cushions and Headgear

|                                    |         |
|------------------------------------|---------|
| EasyFit SilkGel Nasal Mask, Small  | DV97318 |
| EasyFit SilkGel Nasal Mask, Medium | DV97328 |
| EasyFit SilkGel Nasal Mask, Large  | DV97338 |

#### Replacement Cushions

|  |         |
|--|---------|
| EasyFit SilkGel Nasal Mask Cushion, Small  | DV97352 |
| EasyFit SilkGel Nasal Mask Cushion, Medium | DV97353 |
| EasyFit SilkGel Nasal Mask Cushion, Large  | DV97354 |

#### Parts and Accessories

|   |         |
|---|---------|
| EasyFit SilkGel forehead cushion          | DV97270 |
| EasyFit SilkGel forehead support          | DV15856 |
| Ball-and-socket joint set, Nasal          | DV15858 |
| EasyFit SilkGel headgear                  | DV97350 |
| EasyFit SilkGel headgear clips (4 pieces) | DV97305 |
| EasyFit SilkGel rotating sleeve           | DV97288 |
| EasyFit SilkGel oxygen port cap           | DV97308 |
| EasyFit SilkGel tube anchoring strap      | DV97274 |



DeVilbiss Healthcare LLC | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)